




PCYXL rabbit pAb


Catalog No	YP-Ab-11624
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB;ELISA;IHC
Gene Name	PCYOX1L PSEC0105
Immunogen	Synthesized peptide derived from human PCYXL AA range: 406-456
Specificity	This antibody detects endogenous levels of PCYXL at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000;IHC-p 1:50-300; ELISA 2000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PCYOX1L;prenylcysteine oxidase 1 like;Prenylcysteine oxidase like;PCYXL; PCYOX1L PSEC0105
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted.
Tissue Specificity	
Function	cofactor:FAD.,function:Probable oxidoreductase.,similarity:Belongs to the prenylcysteine oxidase family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night



Immunohistochemical analysis of paraffin-embedded human Breast cancer. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA, pH9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200(room temperature, 45min).