




CJ113 rabbit pAb

Catalog No	YP-Ab-11641
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	C10orf113
Immunogen	Synthesized peptide derived from human CJ113 AA range: 73-123
Specificity	This antibody detects endogenous levels of CJ113 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Nebulette;LNEBL;C10orf113;Actin-binding Z-disk protein;Actin binding Z disk protein;CJ113
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	caution:Encoded in intron of the NEBL gene. Confirmed by only 1 EST (AI806552).,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HEK293 cells,
primary antibody was diluted at 1:1000, 4° over night