




T131L rabbit pAb

Catalog No	YP-Ab-11650
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	KIAA0922 TMEM131L
Immunogen	Synthesized peptide derived from human T131L AA range: 321-371
Specificity	This antibody detects endogenous levels of T131L at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KIAA0922;Transmembrane protein 131 like;Transmembrane protein 131-like; T131L;KIAA0922 TMEM131L
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Single-pass type I membrane protein . Cytoplasm . During intrathymic development, resides in punctate cytoplasmic structures in DN1 and DN2 cells. In DN3 cells, found in large crescent-shaped membrane structures, which preferentially localize in cell-to-cell contact zones. .; [Isoform 1]: Endoplasmic reticulum . Transmembrane localization is essential for Wnt signaling inhibition. .; [Isoform 5]: Cytoplasm . Scattered throughout the cytoplasm in small-sized punctate structures. .
Tissue Specificity	Expressed in thymocytes.
Function	similarity:Belongs to the TMEM131 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from HCT116 cells, primary antibody was diluted at 1:1000, 4° over night