




DJC11 rabbit pAb

Catalog No	YP-Ab-11654
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	DNAJC11
Immunogen	Synthesized peptide derived from human DJC11 AA range: 206-256
Specificity	This antibody detects endogenous levels of DJC11 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DNAJC11;DJC11
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Mitochondrion . Isoforms show differential submitochondrial localization. A 57 kDa form (potentially isoform 3) shows either mitochondrial matrix or innermembrane (IM) localization, possibly anchored to the IM facing the matrix. A 35 kDa form behaved either as an inner membrane space (IMS) or an IM protein exposed to the IMS. .; [Isoform 1]: Mitochondrion outer membrane ; Peripheral membrane protein .
Tissue Specificity	
Function	similarity:Belongs to the DNAJC11 family.,similarity:Contains 1 J domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from HeLa cells,
primary antibody was diluted at 1:1000, 4°over night