




F210A rabbit pAb

Catalog No	YP-Ab-11655
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	FAM210A C18orf19
Immunogen	Synthesized peptide derived from human F210A AA range: 176-226
Specificity	This antibody detects endogenous levels of F210A at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C9orf10; Constitutive coactivator of PPAR-gamma-like protein 1; KIAA0183; OSSA; Oxidative stress-; F210A; FAM210A C18orf19
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Single-pass membrane protein . Mitochondrion . Cytoplasm .
Tissue Specificity	
Function	similarity: Contains 1 DUF1279 domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from 3T3 cells, primary antibody was diluted at 1:1000, 4° over night