




R3HD1 rabbit pAb

Catalog No	YP-Ab-11658
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	R3HDM1 KIAA0029 R3HDM
Immunogen	Synthesized peptide derived from human R3HD1 AA range: 327-377
Specificity	This antibody detects endogenous levels of R3HD1 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	R3HD1;R3HDM1 KIAA0029 R3HDM
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	similarity:Contains 1 R3H domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HEK293 cells,
primary antibody was diluted at 1:1000, 4° over night