




CCD41 rabbit pAb

Catalog No	YP-Ab-11662
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	CCDC41
Immunogen	Synthesized peptide derived from human CCD41 AA range: 321-371
Specificity	This antibody detects endogenous levels of CCD41 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CCDC41; CEP83; NY REN 58; CCD41
Observed Band	
Calculated Molecular Weight	76kD
Cell Pathway	Cytoplasm, cytoskeleton, microtubule organizing center, centrosome, centriole . Localizes specifically to the distal appendage region of the centriole, which anchors the mother centriole to the plasma membrane. Localizes to centrioles at all stages of the cell cycle, including mitosis.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HpeG2 cells,
primary antibody was diluted at 1:1000, 4°over night