



PYRD1 rabbit pAb

Catalog No	YP-Ab-11691
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB; ELISA; IHC
Gene Name	PYROXD1
Immunogen	Synthesized peptide derived from human PYRD1 AA range: 374-424
Specificity	This antibody detects endogenous levels of PYRD1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000; IHC-p 1:50-300; ELISA 2000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	pyridine nucleotide disulphide oxidoreductase domain 1; PYRD1; PYROXD1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus . Cytoplasm . Cytoplasm, myofibril, sarcomere .
Tissue Specificity	
Function	cofactor: Binds 1 FAD per subunit., similarity: Belongs to the class-I pyridine nucleotide-disulfide oxidoreductase family. PYROXD1 subfamily.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

Western blot analysis of lysates from HpeG2 cells, primary antibody was diluted at 1:1000, 4° over night

Immunohistochemical analysis of paraffin-embedded human Gastric adenocarcinoma. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA, pH9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200(room temperature, 45min).