




TRI44 rabbit pAb

Catalog No	YP-Ab-11694
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	TRIM44 DIPB
Immunogen	Synthesized peptide derived from human TRI44 AA range: 42-92
Specificity	This antibody detects endogenous levels of TRI44 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DIPB; HSA249128; MC7; Protein DIPB; TRIM44; tripartite motif containing 44; TRI44; TRIM44 DIPB
Observed Band	
Calculated Molecular Weight	
Cell Pathway	intracellular,
Tissue Specificity	Highly expressed in testis.
Function	function: May play a role in the process of differentiation and maturation of neuronal cells., similarity: Contains 1 B box-type zinc finger.,
Background	This gene encodes a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, namely a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from Hela cells,
primary antibody was diluted at 1:1000, 4° over night