




RUSD1 rabbit pAb

Catalog No	YP-Ab-11707
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	RPUUSD1 C16orf40 RLUC
Immunogen	Synthesized peptide derived from human RUSD1 AA range: 58-108
Specificity	This antibody detects endogenous levels of RUSD1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RLUC;C16orf40;RUSD1;RPUUSD1 C16orf40 RLUC
Observed Band	
Calculated Molecular Weight	
Cell Pathway	mitochondrion,
Tissue Specificity	
Function	similarity:Belongs to the pseudouridine synthase rluA family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from 3T3 cells, primary antibody was diluted at 1:1000, 4° over night