




TMM79 rabbit pAb

Catalog No	YP-Ab-11715
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	TMEM79
Immunogen	Synthesized peptide derived from human C-terminal TMM79
Specificity	This antibody detects endogenous levels of TMM79 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TMM79; TMEM79
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Lysosome . Golgi apparatus, trans-Golgi network . Membrane ; Multi-pass membrane protein . Colocalized with TGOLN2 in the trans-Golgi network. Colocalized with LAMP1 in the lysosome (By similarity). .
Tissue Specificity	Expressed in the epidermis of the skin. Expressed in epithelial cells of the outermost layer of the stratum granulosum (SG) and hair follicles (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MCF-7 cells,
primary antibody was diluted at 1:1000, 4°over night