




TDRD9 rabbit pAb

Catalog No	YP-Ab-11733
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	TDRD9 C14orf75
Immunogen	Synthesized peptide derived from human TDRD9 AA range: 335-385
Specificity	This antibody detects endogenous levels of TDRD9 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TDRD9; TDRD9 C14orf75
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Nucleus . Component of the nuage, also named P granule, a germ-cell-specific organelle required to repress transposon activity during meiosis. Specifically localizes to piP-bodies, a subset of the nuage which contains secondary piRNAs. PIWIL2 is required for its localization to piP-bodies. .
Tissue Specificity	
Function	similarity: Contains 1 helicase ATP-binding domain., similarity: Contains 1 helicase C-terminal domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HpeG2 cells,
primary antibody was diluted at 1:1000, 4°over night