



GAGB1 rabbit pAb


Catalog No	YP-Ab-11735
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	PAGE1 GAGE9 GAGEB1
Immunogen	Synthesized peptide derived from human GAGB1 AA range: 28-78
Specificity	This antibody detects endogenous levels of GAGB1 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GAGB1;PAGE1 GAGE9 GAGEB1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Isolated from prostate cancer cell lines; expression associated with progression to androgen insensitive phenotype. Expressed in normal testis and at lower level in normal placenta.
Function	function:Unknown.,similarity:Belongs to the GAGE family.,tissue specificity:Isolated from prostate cancer cell lines; expression associated with progression to androgen insensitive phenotype. Expressed in normal testis and at lower level in normal placenta.,
Background	This gene belongs to a family of genes that are expressed in a variety of tumors but not in normal tissues, except for the testis. Unlike the other gene family members, this gene does not encode an antigenic peptide. Nothing is presently known about the function of this protein. [provided by RefSeq, Jul 2008],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from CACO2 cells, primary antibody was diluted at 1:1000, 4°over night