



FAT3 rabbit pAb

Catalog No	YP-Ab-11741
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	IHC; IF
Gene Name	FAT3 CDHF15 KIAA1989
Immunogen	Synthesized peptide derived from human FAT3 AA range: 896-946
Specificity	This antibody detects endogenous levels of FAT3 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:50-200. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FAT3; FAT3 CDHF15 KIAA1989
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Single-pass type I membrane protein .
Tissue Specificity	Expressed in ES cells, primitive neuroectoderm, fetal brain, infant brain, adult neural tissues and prostate.
Function	function: May play a role in the interactions between neurites derived from specific subsets of neurons during development., similarity: Contains 1 laminin G-like domain., similarity: Contains 33 cadherin domains., similarity: Contains 4 EGF-like domains., tissue specificity: Expressed in ES cells, primitive neuroectoderm, fetal brain, infant brain, adult neural tissues and prostate.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images