



LACB2 rabbit pAb

Catalog No	YP-Ab-11750
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB; ELISA; IHC
Gene Name	LACTB2 CGI-83
Immunogen	Synthesized peptide derived from human LACB2 AA range: 19-69
Specificity	This antibody detects endogenous levels of LACB2 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000; IHC-p 1:50-300; ELISA 2000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EC:3.1.27.-; CGI-83; CGI83; CGI 83; Beta-lactamase-like protein 2; LACB2; LACTB2 CGI-83
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Mitochondrion matrix .
Tissue Specificity	
Function	similarity: Belongs to the metallo-beta-lactamase superfamily. Glyoxalase II family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

Western blot analysis of lysates from 293T cells, primary antibody was diluted at 1:1000, 4° over night

Immunohistochemical analysis of paraffin-embedded human liver cancer. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA, pH9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200(room temperature, 45min).