



LEAP2 rabbit pAb

Catalog No	YP-Ab-11756
Isotype	IgG
Reactivity	Human; Mouse
Applications	IHC;IF
Gene Name	LEAP2
Immunogen	Synthesized peptide derived from human LEAP2 AA range: 21-71
Specificity	This antibody detects endogenous levels of LEAP2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:50-200. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LEAP2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted.
Tissue Specificity	
Function	function:Has an antimicrobial activity.,similarity:Belongs to the LEAP2 family.,
Background	This gene encodes a cysteine-rich cationic antimicrobial peptide that is expressed predominantly in the liver. The mature peptide has activity against gram-positive bacteria and yeasts. [provided by RefSeq, Sep 2014],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 30min).