



K1551 rabbit pAb

Catalog No	YP-Ab-11760
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB;IHC
Gene Name	KIAA1551 C12orf35
Immunogen	Synthesized peptide derived from human K1551 AA range: 262-312
Specificity	This antibody detects endogenous levels of K1551 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000;IHC-p 1:50-300
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C12orf35;KIAA1551;Retroelement silencing factor 1;K1551;KIAA1551 C12orf35
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus . Localizes around gamma-tubulin during M phase. .
Tissue Specificity	
Function	sequence caution:Contaminating sequence.,sequence caution:Contaminating sequence. Potential poly-A sequence.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

Western blot analysis of lysates from VEC cells, primary antibody was diluted at 1:1000, 4° over night

Immunohistochemical analysis of paraffin-embedded human liver cancer. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA, pH9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200(room temperature, 45min).