




P33MX rabbit pAb

Catalog No	YP-Ab-11761
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	KIAA1191 P33MONOX
Immunogen	Synthesized peptide derived from human P33MX AA range: 186-236
Specificity	This antibody detects endogenous levels of P33MX at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Putative monooxygenase p33MONOX; P33MONOX; Flavin monooxygenase motif-containing protein of 33 kDa; EC:1.-.-.-; P33MX; KIAA1191 P33MONOX
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm .
Tissue Specificity	Down-regulated in the occipital lobe of an early stage Alzheimer disease patients.
Function	similarity: Belongs to the UPF0498 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from SH-SY5Y cells,
primary antibody was diluted at 1:1000, 4° over night