




HFM1 rabbit pAb

Catalog No	YP-Ab-11763
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	HFM1 SEC3D1
Immunogen	Synthesized peptide derived from human HFM1 AA range: 943-993
Specificity	This antibody detects endogenous levels of HFM1 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SEC3D1;HFM1;HFM1 SEC3D1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Preferentially expressed in testis and ovary.
Function	similarity:Belongs to the helicase family. SKI2 subfamily.,similarity:Contains 1 helicase ATP-binding domain.,similarity:Contains 1 helicase C-terminal domain.,similarity:Contains 1 SEC63 domain.,
Background	The protein encoded by this gene is thought to be an ATP-dependent DNA helicase and is expressed mainly in germ-line cells. Defects in this gene are a cause of premature ovarian failure 9 (POF9). [provided by RefSeq, Apr 2014],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from U2OS cells,
primary antibody was diluted at 1:1000, 4°over night