




CTGE3 rabbit pAb

Catalog No	YP-Ab-11769
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	CTAGE3P CTAGE3
Immunogen	Synthesized peptide derived from human CTGE3 AA range: 76-126
Specificity	This antibody detects endogenous levels of CTGE3 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CTGE3;CTAGE3P CTAGE3
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Expressed in normal tissues including colon, mammary gland, ovary, placenta, stomach and testis, as well as several fetal tissues.
Function	disease:Tumor-associated antigen found in several cutaneous T-cell lymphoma (CTCL). Expressed in 4 of 11 Mycosis fungoides tissues and in a T-zone lymphoma.,function:Tumor-associated antigen.,similarity:Belongs to the cTAGE family.,tissue specificity:Expressed in normal tissues including colon, mammary gland, ovary, placenta, stomach and testis, as well as several fetal tissues.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from A549 cells,
primary antibody was diluted at 1:1000, 4° over night