




F220A rabbit pAb

Catalog No	YP-Ab-11783
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	FAM220A C7orf70 SIPAR
Immunogen	Synthesized peptide derived from human F220A AA range: 186-236
Specificity	This antibody detects endogenous levels of F220A at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	F220A;FAM220A C7orf70 SIPAR
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	
Function	function:May negatively regulate STAT3.,subunit:Interacts with STAT3.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from PC-12 cells,
primary antibody was diluted at 1:1000, 4° over night