




GOG8H rabbit pAb

Catalog No	YP-Ab-11787
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	GOLGA8H
Immunogen	Synthesized peptide derived from human GOG8H AA range: 478-528
Specificity	This antibody detects endogenous levels of GOG8H at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GOG8H;GOLGA8H
Observed Band	
Calculated Molecular Weight	
Cell Pathway	atus,
Tissue Specificity	
Function	caution:The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data., similarity:Belongs to the GOLGA8 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from A549 cells,
primary antibody was diluted at 1:1000, 4° over night