




PAMR1 rabbit pAb

Catalog No	YP-Ab-11797
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	PAMR1 RAMP FP938 UNQ699/PRO1344
Immunogen	Synthesized peptide derived from human PAMR1 AA range: 57-107
Specificity	This antibody detects endogenous levels of PAMR1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DKFZP586H2123;FP938;Inactive serine protease PAMR1;PAMR1;RAMP; UNQ699/PRO1344;PAMR1 RAMP FP938 UNQ699/PRO1344
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HpeG2 cells,
primary antibody was diluted at 1:1000, 4° over night