




TB22A rabbit pAb

Catalog No	YP-Ab-11799
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	TBC1D22A C22orf4
Immunogen	Synthesized peptide derived from human TB22A AA range: 49-99
Specificity	This antibody detects endogenous levels of TB22A at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TBC1 domain family;member 22A;TBC1 domain family member 22A; HSC79E021;C22orf4;TB22A;TBC1D22A C22orf4
Observed Band	
Calculated Molecular Weight	
Cell Pathway	intracellular,endomembrane system,
Tissue Specificity	
Function	function:May act as a GTPase-activating protein for Rab family protein(s)., similarity:Contains 1 Rab-GAP TBC domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from PC-12 cells,
primary antibody was diluted at 1:1000, 4° over night