




RNF10 rabbit pAb


Catalog No	YP-Ab-11807
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB; ELISA; IHC
Gene Name	RNF10 KIAA0262 RIE2
Immunogen	Synthesized peptide derived from human RNF10 AA range: 294-344
Specificity	This antibody detects endogenous levels of RNF10 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000; IHC-p 1:50-300; ELISA 2000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RING finger protein 10; RIE2; KIAA0262; EC:2.3.2.27; E3 ubiquitin-protein ligase RNF10; RNF10; RNF10 KIAA0262 RIE2
Observed Band	89kD
Calculated Molecular Weight	89kD
Cell Pathway	Cytoplasm . Nucleus .
Tissue Specificity	
Function	similarity: Belongs to the RNF10 family., similarity: Contains 1 RING-type zinc finger., subunit: Interacts with MEOX2.,
Background	The protein encoded by this gene contains a ring finger motif, which is known to be involved in protein-protein interactions. The specific function of this protein has not yet been determined. EST data suggests the existence of multiple alternatively spliced transcript variants, however, their full length nature is not known. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from Jarkat cells, primary antibody was diluted at 1:1000, 4° over night



Immunohistochemical analysis of paraffin-embedded human Squamous cell carcinoma of lung. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA, pH9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200(room temperature, 45min).