




CA170 rabbit pAb

Catalog No	YP-Ab-11826
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	C1orf170
Immunogen	Synthesized peptide derived from human CA170 AA range: 434-484
Specificity	This antibody detects endogenous levels of CA170 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MYG;Gamm1;MST024;MSTP024;C12orf10;CA170;C1orf170
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Nucleus . Shows a nuclear localization in the presence of PPARGC1A. .
Tissue Specificity	Muscle-specific expression is increased by endurance exercise.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MDA-MB cells,
primary antibody was diluted at 1:1000, 4° over night