



CEL3B rabbit pAb

Catalog No	YP-Ab-11833
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	CELA3B ELA3B
Immunogen	Synthesized peptide derived from human CEL3B AA range: 206-256
Specificity	This antibody detects endogenous levels of CEL3B at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Protease E;Elastase-3B;Elastase IIIB;EC:3.4.21.70;Chymotrypsin-like elastase family member 3B;CEL3B;CELA3B ELA3B
Observed Band	
Calculated Molecular Weight	
Cell Pathway	extracellular space,
Tissue Specificity	Pancreas. Not detectable in keratinocytes.
Function	
Background	Elastases form a subfamily of serine proteases that hydrolyze many proteins in addition to elastin. Humans have six elastase genes which encode the structurally similar proteins elastase 1, 2, 2A, 2B, 3A, and 3B. Unlike other elastases, elastase 3B has little elastolytic activity. Like most of the human elastases, elastase 3B is secreted from the pancreas as a zymogen and, like other serine proteases such as trypsin, chymotrypsin and kallikrein, it has a digestive function in the intestine. Elastase 3B preferentially cleaves proteins after alanine residues. Elastase 3B may also function in the intestinal transport and metabolism of cholesterol. Both elastase 3A and elastase 3B have been



referred to as protease E and as elastase 1, and excretion of this protein in fecal material is frequently used as a measure of pancreatic function in clinical assays. [provided by RefSeq, May 2009],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from SH-SY5Y cells, primary antibody was diluted at 1:1000, 4° over night