




GNT2A rabbit pAb

Catalog No	YP-Ab-11834
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	GCNT2 GCNT5 II NACGT1
Immunogen	Synthesized peptide derived from human GNT2A AA range: 191-241
Specificity	This antibody detects endogenous levels of GNT2A at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CCAT;EC:2.4.1.150;GCNT2C;GCNT5;l-branching enzyme;GNT2A;GCNT2 GCNT5 II NACGT1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Golgi apparatus membrane ; Single-pass type II membrane protein .
Tissue Specificity	[Isoform B]: Expressed in lens epithelium cells. ; [Isoform C]: Expressed in reticulocytes.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HeLa cells,
primary antibody was diluted at 1:1000, 4° over night