




VWA8 rabbit pAb

Catalog No	YP-Ab-11839
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	VWA8 KIAA0564
Immunogen	Synthesized peptide derived from human VWA8 AA range: 1853-1903
Specificity	This antibody detects endogenous levels of VWA8 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RP97;P7BP2;KIAA0564;VWA8;VWA8 KIAA0564
Observed Band	
Calculated Molecular Weight	
Cell Pathway	[Isoform 1]: Mitochondrion .; [Isoform 2]: Mitochondrion . Peroxisome . Localizes to peroxisomes in a PEX7-dependent manner. .
Tissue Specificity	
Function	similarity:Contains 1 VWFA domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from A431 cells, primary antibody was diluted at 1:1000, 4° over night