




TM192 rabbit pAb


Catalog No	YP-Ab-11844
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB; IHC
Gene Name	TMEM192
Immunogen	Synthesized peptide derived from human TM192 AA range: 27-77
Specificity	This antibody detects endogenous levels of TM192 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000; IHC-p 1:50-300
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FLJ38482; TMEM192; transmembrane protein 192; TM192
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Lysosome membrane ; Multi-pass membrane protein . Late endosome .
Tissue Specificity	Strongly expressed in kidney, liver, lung and pancreas.
Function	PTM: Phosphorylated upon DNA damage, probably by ATM or ATR., similarity: Belongs to the TMEM192 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MCF-7 cells, primary antibody was diluted at 1:1000, 4° over night



Immunohistochemical analysis of paraffin-embedded human liver cancer. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA, pH9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200(room temperature, 45min).