



NEUR2 rabbit pAb

Catalog No	YP-Ab-11852
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	NEU2
Immunogen	Synthesized peptide derived from human NEUR2 AA range: 69-119
Specificity	This antibody detects endogenous levels of NEUR2 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SIAL2; NEU2; NEUR2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, cytosol .
Tissue Specificity	Expressed in skeletal muscle, fetal liver and embryonic carcinoma cell line NT2-D1.
Function	catalytic activity: Hydrolysis of alpha-(2->3)-, alpha-(2->6)-, alpha-(2->8)- glycosidic linkages of terminal sialic acid residues in oligosaccharides, glycoproteins, glycolipids, colominic acid and synthetic substrates., function: Hydrolyzes sialylated compounds., similarity: Belongs to the glycosyl hydrolase 33 family., similarity: Contains 2 BNR repeats., tissue specificity: Expressed in skeletal muscle, fetal liver and embryonic carcinoma cell line NT2D1.,
Background	This gene belongs to a family of glycohydroltic enzymes which remove sialic acid residues from glycoproteins and glycolipids. Expression studies in COS7 cells confirmed that this gene encodes a functional sialidase. Its cytosolic localization was demonstrated by cell fractionation experiments. [provided by

RefSeq, Jul 2008],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from MDA-MB cells, primary antibody was diluted at 1:1000, 4°over night