




TAC2N rabbit pAb

Catalog No	YP-Ab-11855
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	TC2N C14orf47 MTAC2D1
Immunogen	Synthesized peptide derived from human TAC2N AA range: 304-354
Specificity	This antibody detects endogenous levels of TAC2N at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TAC2N;TC2N C14orf47 MTAC2D1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	
Function	similarity:Contains 1 C2 domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from CACO2 cells,
primary antibody was diluted at 1:1000, 4° over night