




PZRN4 rabbit pAb

Catalog No	YP-Ab-11859
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	PDZRN4 LNX4 SEMCAP3L
Immunogen	Synthesized peptide derived from human PZRN4 AA range: 287-337
Specificity	This antibody detects endogenous levels of PZRN4 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PZRN4;PDZRN4 LNX4 SEMCAP3L
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	similarity:Contains 1 RING-type zinc finger.,similarity:Contains 2 PDZ (DHR) domains.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from A549 cells,
primary antibody was diluted at 1:1000, 4° over night