




CL076 rabbit pAb

Catalog No	YP-Ab-11863
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	C12orf76
Immunogen	Synthesized peptide derived from human CL076 AA range: 59-109
Specificity	This antibody detects endogenous levels of CL076 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CL076;C12orf76
Observed Band	
Calculated Molecular Weight	
Cell Pathway	integral component of plasma membrane,integral component of membrane,
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MCF-7 cells,
primary antibody was diluted at 1:1000, 4° over night