




CU082 rabbit pAb

Catalog No	YP-Ab-11871
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	IHC;IF
Gene Name	LINC00310 C21orf82 NCRNA00310
Immunogen	Synthesized peptide derived from human CU082 AA range: 12-62
Specificity	This antibody detects endogenous levels of CU082 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:50-200. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CU082;LINC00310 C21orf82 NCRNA00310
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Widely expressed; not found in breast.
Function	tissue specificity:Widely expressed; not found in breast.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 30min).