




VNN3 rabbit pAb

Catalog No	YP-Ab-11873
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	VNN3
Immunogen	Synthesized peptide derived from human VNN3 AA range: 87-137
Specificity	This antibody detects endogenous levels of VNN3 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	VNN3
Observed Band	
Calculated Molecular Weight	
Cell Pathway	extracellular space,plasma membrane,anchored component of membrane,
Tissue Specificity	
Function	function:Probable hydrolase.,similarity:Belongs to the CN hydrolase family. BTD/ VNN subfamily.,similarity:Contains 1 CN hydrolase domain.,
Background	This gene is the central gene in a cluster of three vanin genes on chromosome 6q23-q24. Extensive alternative splicing has been described; the two most common variants are represented as RefSeqs. [provided by RefSeq, Apr 2014],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HuvEc cells,
primary antibody was diluted at 1:1000, 4° over night