




CB083 rabbit pAb

Catalog No	YP-Ab-11875
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	C2orf83
Immunogen	Synthesized peptide derived from human CB083 AA range: 56-106
Specificity	This antibody detects endogenous levels of CB083 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CB083;C2orf83
Observed Band	
Calculated Molecular Weight	
Cell Pathway	integral component of plasma membrane,
Tissue Specificity	
Function	caution:Although related to reduced folate carrier family, it lacks transmembrane domains, suggesting that it has probably no transporter function.,sequence caution:Translated as Trp.,similarity:Belongs to the reduced folate carrier (RFC) transporter (TC 2.A.48) family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from U2OS cells,
primary antibody was diluted at 1:1000, 4° over night