




VGLL1 rabbit pAb

Catalog No	YP-Ab-11877
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	VGLL1 TDU
Immunogen	Synthesized peptide derived from human VGLL1 AA range: 84-134
Specificity	This antibody detects endogenous levels of VGLL1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Protein TONDU;TDU;vestigial like 1 (Drosophila);Vgl 1;VGL1;VGLL1;VGLL1 TDU
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus.
Tissue Specificity	
Function	function:May act as a specific coactivator for the mammalian TEFs., similarity:Belongs to the vestigial family.,subunit:Interacts with TEFs.,
Background	The protein encoded by this gene binds proteins of the TEA domain family of transcription factors (TEFs) through the Vg (vestigial) homology region found in its N-terminus. It may thus function as a specific coactivator for the mammalian TEFs. [provided by RefSeq, Sep 2009],
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from Hela cells, primary antibody was diluted at 1:1000, 4° over night