




DEN2C rabbit pAb

Catalog No	YP-Ab-11880
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	DENND2C
Immunogen	Synthesized peptide derived from human DEN2C AA range: 763-813
Specificity	This antibody detects endogenous levels of DEN2C at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DEN2C;DENND2C
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	similarity:Contains 1 dDENN domain.,similarity:Contains 1 DENN domain., similarity:Contains 1 uDENN domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from Hela cells,
primary antibody was diluted at 1:1000, 4° over night