




PNM6A rabbit pAb

Catalog No	YP-Ab-11881
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	PNMA6A
Immunogen	Synthesized peptide derived from human PNM6A AA range: 171-221
Specificity	This antibody detects endogenous levels of PNM6A at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PNM6A;PNMA6A
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Expressed in the brain.
Function	similarity:Belongs to the PNMA family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from THP-1 cells,
primary antibody was diluted at 1:1000, 4° over night