




NSUN7 rabbit pAb

Catalog No	YP-Ab-11882
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	NSUN7
Immunogen	Synthesized peptide derived from human NSUN7 AA range: 139-189
Specificity	This antibody detects endogenous levels of NSUN7 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FLJ14001;NSUN7
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	function:May have S-adenosyl-L-methionine-dependent methyl-transferase activity .,similarity:Belongs to the methyltransferase superfamily. RsmB/NOP family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MDA-MB cells,
primary antibody was diluted at 1:1000, 4° over night