




CG013 rabbit pAb

Catalog No	YP-Ab-11888
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	C7orf13
Immunogen	Synthesized peptide derived from human CG013 AA range: 133-183
Specificity	This antibody detects endogenous levels of CG013 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CG013;C7orf13
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Membrane ; Single-pass membrane protein .
Tissue Specificity	Expressed only in testis.
Function	tissue specificity:Expressed only in testis.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MCF-7 cells,
primary antibody was diluted at 1:1000, 4° over night