




CC124 rabbit pAb

Catalog No	YP-Ab-11907
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	CCDC124
Immunogen	Synthesized peptide derived from human CC124 AA range: 56-106
Specificity	This antibody detects endogenous levels of CC124 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	coiled coil domain containing 124;Coiled-coil domain-containing protein 124; CC124;CCDC124
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, cytoskeleton, microtubule organizing center, centrosome . Midbody . Colocalizes with gamma-tubulin at interphase, prophase, metaphase, and anaphase. Relocates from centrosome to midbody at telophase.
Tissue Specificity	Ubiquitously expressed.
Function	similarity:Belongs to the CCDC124 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from A431 cells,
primary antibody was diluted at 1:1000, 4° over night