




HMN2 rabbit pAb

Catalog No	YP-Ab-11908
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	IHC;IF
Gene Name	MTRNR2L2
Immunogen	Synthesized peptide derived from human HMN2 AA range: 1-24
Specificity	This antibody detects endogenous levels of HMN2 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:50-200. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	H11;HMN2;CMT2L;DHMN2;E2IG1;HMN2A;HSP22;MTRNR2L2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted . Cytoplasm .
Tissue Specificity	Highly expressed in testis. Also expressed in kidney, heart, skeletal muscles and brain.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 30min).