




TIDC1 rabbit pAb


Catalog No	YP-Ab-11913
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB; ELISA; IHC
Gene Name	TIMMDC1 C3orf1 UNQ247/PRO284
Immunogen	Synthesized peptide derived from human TIDC1 AA range: 74-124
Specificity	This antibody detects endogenous levels of TIDC1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000; IHC-p 1:50-300; ELISA 2000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C3orf1; Protein M5 14; TIMMDC1; Transmembrane protein C3orf1; TIDC1; TIMMDC1 C3orf1 UNQ247/PRO284
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Mitochondrion membrane ; Multi-pass membrane protein .
Tissue Specificity	Generalized expression enhanced in heart and skeletal muscle.
Function	tissue specificity: Generalized expression enhanced in heart and skeletal muscle.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from VEC cells, primary antibody was diluted at 1:1000, 4° over night



Immunohistochemical analysis of paraffin-embedded human liver cancer. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA, pH9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200(room temperature, 45min).