




ZCHC7 rabbit pAb

Catalog No	YP-Ab-11915
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	ZCCHC7 HSPC086
Immunogen	Synthesized peptide derived from human ZCHC7 AA range: 4-54
Specificity	This antibody detects endogenous levels of ZCHC7 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AIR1; HSPC086; RP11 397D12.1; ZCCHC7; ZCHC7; ZCCHC7 HSPC086
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus, nucleolus .
Tissue Specificity	
Function	caution: It is uncertain whether Met-1 or Met-2 is the initiator., sequence caution: Contaminating sequence. Potential poly-A sequence., similarity: Contains 4 CCHC-type zinc fingers.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from HeLa cells,
primary antibody was diluted at 1:1000, 4° over night