




MBLC1 rabbit pAb


Catalog No	YP-Ab-11920
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB; IHC
Gene Name	MBLAC1
Immunogen	Synthesized peptide derived from human MBLC1 AA range: 1-51
Specificity	This antibody detects endogenous levels of MBLC1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000; IHC-p 1:50-300
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MBLC1; MBLAC1
Observed Band	29kD
Calculated Molecular Weight	29kD
Cell Pathway	Cytoplasm, cytosol . Nucleus . Localizes in the nucleus during early S-phase of the cell cycle. .
Tissue Specificity	
Function	cofactor: Binds 2 zinc ions per subunit., similarity: Belongs to the metallo-beta-lactamase superfamily. Glyoxalase II family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MCF-7 cells, primary antibody was diluted at 1:1000, 4° over night



Immunohistochemical analysis of paraffin-embedded human meningioma. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA, pH9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200(room temperature, 45min).