




ERI2 rabbit pAb

Catalog No	YP-Ab-11928
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	ERI2 EXOD1 KIAA1504
Immunogen	Synthesized peptide derived from human ERI2 AA range: 430-480
Specificity	This antibody detects endogenous levels of ERI2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ERI2;ERI2 EXOD1 KIAA1504
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from MDA-MB cells,
primary antibody was diluted at 1:1000, 4° over night