




OVOS2 rabbit pAb

Catalog No	YP-Ab-11930
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	OVOS2
Immunogen	Synthesized peptide derived from human OVOS2 AA range: 660-710
Specificity	This antibody detects endogenous levels of OVOS2 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Ovostatin homolog 2;OVOS2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	
Function	function:Is able to inhibit all four classes of proteinases by a unique 'trapping' mechanism.,similarity:Belongs to the protease inhibitor I39 (alpha-2-macroglobulin) family.,subunit:Homotetramer.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MDA-MB cells,
primary antibody was diluted at 1:1000, 4° over night