




RTP4 rabbit pAb

Catalog No	YP-Ab-11933
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	RTP4 IFRG28
Immunogen	Synthesized peptide derived from human RTP4 AA range: 189-239
Specificity	This antibody detects endogenous levels of RTP4 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RTP4;RTP4 IFRG28
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Single-pass type III membrane protein .
Tissue Specificity	Expressed in circumvallate papillae and testis.
Function	induction:By interferons.,similarity:Belongs to the TMEM7 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from A431 cells,
primary antibody was diluted at 1:1000, 4° over night